

BCL2/adenovirus E1B 19 kDa protein-interacting protein 3

Catalog No: #AP81102

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #AP81102-1 50ug #AP81102-2 100ug #AP81102-3 1mg

Description

Product Name	BCL2/adenovirus E1B 19 kDa protein-interacting protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Nip3
Accession No.	Uniprot ID:O55003
Uniprot	O55003
GeneID	12176;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Bnip3

Note: This product is for in vitro research use only