

Apoptosis regulator Bcl-2

Catalog No: #AP81103



Package Size: #AP81103-1 50ug #AP81103-2 100ug #AP81103-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Apoptosis regulator Bcl-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Bcl-2
Accession No.	Uniprot ID:P10417
Uniprot	P10417
GeneID	12043;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Bcl2

Note: This product is for in vitro research use only