

B-cell linker protein

Catalog No: #AP81105



Package Size: #AP81105-1 50ug #AP81105-2 100ug #AP81105-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	B-cell linker protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	B-cell adapter containing a SH2 domain protein,B-cell adapter containing a Src homology 2 domain protein,Cytoplasmic adapter protein,Src homology 2 domain-containing leukocyte protein of 65 kDa,SLP-65,BASH,SLP65
Accession No.	Uniprot ID:Q8WV28
Uniprot	Q8WV28
GeneID	29760;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BLNK

Note: This product is for in vitro research use only