

# Luciferase Mouse Monoclonal Antibody(6B8)

Catalog No: #44173

Package Size: #44173-1 50ul #44173-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Luciferase Mouse Monoclonal Antibody(6B8)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	6B8
Isotype	IgG1
Applications	WB
Species Reactivity	Dm
Specificity	Luciferase Mouse Monoclonal Antibody(6B8) detects endogenous levels of Luciferase
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	Luciferase
Other Names	Luciferase
Accession No.	Swiss-Prot:P08659Gene ID:
Uniprot	P08659
Calculated MW	60kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western blotting: 1:1000-1:2000

## Images



Western Blot analysis of Luciferase protein using antibody diluted at 1:1000

---

Note: This product is for in vitro research use only