

Collagen type IV alpha-3-binding protein

Catalog No: #AP81111



Package Size: #AP81111-1 50ug #AP81111-2 100ug #AP81111-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Collagen type IV alpha-3-binding protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Ceramide transfer protein,Goodpasture antigen-binding protein,GPBP,START domain-containing protein 11,StARD11,StAR-related lipid transfer protein 11,Stard11
Accession No.	Uniprot ID:Q9EQG9
Uniprot	Q9EQG9
GeneID	68018;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Col4a3bp

Note: This product is for in vitro research use only