

# Brain-derived neurotrophic factor

Catalog No: #AP81115



Package Size: #AP81115-1 50ug #AP81115-2 100ug #AP81115-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Brain-derived neurotrophic factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	BDNF
Accession No.	Uniprot ID:Q7YRB4
Uniprot	Q7YRB4
GeneID	403461;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: BDNF

Note: This product is for in vitro research use only