

A disintegrin and metalloproteinase with thrombospondin motifs 4

Catalog No: #AP81120

Package Size: #AP81120-1 50ug #AP81120-2 100ug #AP81120-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	A disintegrin and metalloproteinase with thrombospondin motifs 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ADAM-TS 4,ADMP-1,Aggrecanase-1,KIAA0688,UNQ769/PRO1563
Accession No.	Uniprot ID:O75173
Uniprot	O75173
GeneID	9507;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ADAMTS4

Note: This product is for in vitro research use only