

Apoptosis regulator BAX

Catalog No: #AP81123



Package Size: #AP81123-1 50ug #AP81123-2 100ug #AP81123-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Apoptosis regulator BAX
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Bcl-2-like protein 4,Bcl2-L-4,BCL2L4
Accession No.	Uniprot ID:Q07812
Uniprot	Q07812
GeneID	581;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BAX

Note: This product is for in vitro research use only