

Myc box-dependent-interacting protein 1

Catalog No: #AP81124



Package Size: #AP81124-1 50ug #AP81124-2 100ug #AP81124-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myc box-dependent-interacting protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Amphiphysin II,Amphiphysin-like protein,Bridging integrator 1,SH3 domain-containing protein 9,AmphI,Sh3p9
Accession No.	Uniprot ID:O08539
Uniprot	O08539
GeneID	30948;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Bin1

Note: This product is for in vitro research use only