

Collagen type IV alpha-3-binding protein

Catalog No: #AP81125



Package Size: #AP81125-1 50ug #AP81125-2 100ug #AP81125-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Collagen type IV alpha-3-binding protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Ceramide transfer protein, Goodpasture antigen-binding protein, GPBP, START domain-containing protein 11, StARD11, StAR-related lipid transfer protein 11, STARD11
Accession No.	Uniprot ID:Q9GKI7
Uniprot	Q9GKI7
GeneID	286797;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: COL4A3BP

Note: This product is for in vitro research use only