

## C5a anaphylatoxin chemotactic receptor C5L2

Catalog No: #AP81126



Package Size: #AP81126-1 50ug #AP81126-2 100ug #AP81126-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	C5a anaphylatoxin chemotactic receptor C5L2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	G-protein coupled receptor 77,C5I2
Accession No.	Uniprot ID:Q8BW93
Uniprot	Q8BW93
GeneID	319430;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Gpr77

Note: This product is for in vitro research use only