

Cadherin-5

Catalog No: #AP81127



Package Size: #AP81127-1 50ug #AP81127-2 100ug #AP81127-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cadherin-5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Vascular endothelial cadherin,VE-cadherin
Accession No.	Uniprot ID:Q8AYD0
Uniprot	Q8AYD0
GeneID	374068;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CDH5

Note: This product is for in vitro research use only