

Calbindin

Catalog No: #AP81130



Package Size: #AP81130-1 50ug #AP81130-2 100ug #AP81130-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Calbindin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Calbindin D28,D-28K,Spot 35 protein,Vitamin D-dependent calcium-binding protein, avian-type
Accession No.	Uniprot ID:P07171
Uniprot	P07171
GeneID	83839;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Calb1

Note: This product is for in vitro research use only