

Breast cancer type 2 susceptibility protein

Catalog No: #AP81131



Package Size: #AP81131-1 50ug #AP81131-2 100ug #AP81131-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Breast cancer type 2 susceptibility protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Fanconi anemia group D1 protein,FACD,FANCD1
Accession No.	Uniprot ID:P51587
Uniprot	P51587
GeneID	675;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BRCA2

Note: This product is for in vitro research use only