

Tyrosine-protein kinase BTK

Catalog No: #AP81132



Package Size: #AP81132-1 50ug #AP81132-2 100ug #AP81132-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tyrosine-protein kinase BTK
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Agammaglobulinemia tyrosine kinase,ATK,B-cell progenitor kinase,BPK,Bruton tyrosine kinase,Kinase EMB,Bpk
Accession No.	Uniprot ID:P35991
Uniprot	P35991
GeneID	12229;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Btk

Note: This product is for in vitro research use only