

Complement C1q tumor necrosis factor-related protein

3

Catalog No: #AP81133

Package Size: #AP81133-1 50ug #AP81133-2 100ug #AP81133-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Complement C1q tumor necrosis factor-related protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Collagenous repeat-containing sequence 26 kDa protein,CORS26,Secretory protein CORS26,CTRP3,UNQ753/PRO1484
Accession No.	Uniprot ID:Q9BXJ4
Uniprot	Q9BXJ4
GeneID	114899;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:C1QTNF3

Note: This product is for in vitro research use only