

Complement C1q subcomponent subunit C

Catalog No: #AP81140



Package Size: #AP81140-1 50ug #AP81140-2 100ug #AP81140-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Complement C1q subcomponent subunit C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	C1QG
Accession No.	Uniprot ID:P02747
Uniprot	P02747
GeneID	714;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:C1QC

Note: This product is for in vitro research use only