

C-C motif chemokine 4

Catalog No: #AP81142



Package Size: #AP81142-1 50ug #AP81142-2 100ug #AP81142-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C motif chemokine 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Small-inducible cytokine A4, SCY4A
Accession No.	Uniprot ID:Q17QA1
Uniprot	Q17QA1
GeneID	414347;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CCL4

Note: This product is for in vitro research use only