

# G1/S-specific cyclin-D3

Catalog No: #AP81143



Package Size: #AP81143-1 50ug #AP81143-2 100ug #AP81143-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	G1/S-specific cyclin-D3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Accession No.	Uniprot ID:Q3MHH5
Uniprot	Q3MHH5
GeneID	540547;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CCND3

Note: This product is for in vitro research use only