

Beclin-1

Catalog No: #AP81148

Package Size: #AP81148-1 50ug #AP81148-2 100ug #AP81148-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beclin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Coiled-coil myosin-like BCL2-interacting protein,Protein GT197
Accession No.	Uniprot ID:Q91XJ1
Uniprot	Q91XJ1
GeneID	114558;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Becn1

Note: This product is for in vitro research use only