

Cyclin B1 Mouse Monoclonal Antibody(1A5)

Catalog No: #44177

Package Size: #44177-1 50ul #44177-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	Cyclin B1 Mouse Monoclonal Antibody(1A5)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	1A5
Isotype	IgG1
Applications	WB IHC
Species Reactivity	Hu Rt Ms
Specificity	Cyclin B1 Mouse Monoclonal Antibody(1A5) detects endogenous levels of CCNB1
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	Cyclin B1
Other Names	CCNB1
Accession No.	Swiss-Prot:P14635Gene ID:891
Uniprot	P14635
GeneID	891;
Calculated MW	55kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western blotting: 1:1000-1:2000

Immunohistochemistry: 1:100-1:200

Images



Western blot analysis of 1) Hela Cell Lysate, 2) 3T3 Cell Lysate, 3) Rat Brain Tissue Lysate using LC3A Mouse mAb diluted at 1:1000.



Immunohistochemical analysis of paraffin-embedded Human Tonsil Tissue using Cyclin B1 Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only