

B-cell linker protein

Catalog No: #AP81152

Package Size: #AP81152-1 50ug #AP81152-2 100ug #AP81152-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | B-cell linker protein |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | B-cell adapter containing a SH2 domain protein,B-cell adapter containing a Src homology 2 domain protein,Cytoplasmic adapter protein,Lymphocyte antigen 57,Src homology 2 domain-containing leukocyte protein of 65 kDa,Slp-65,Bash,Ly57,Slp65 |
| Accession No. | Uniprot ID:Q9QUN3 |
| Uniprot | Q9QUN3 |
| GeneID | 17060; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Blnk

Note: This product is for in vitro research use only