

Bone morphogenetic protein 6

Catalog No: #AP81153



Package Size: #AP81153-1 50ug #AP81153-2 100ug #AP81153-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Bone morphogenetic protein 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	BMP-6,VG-1-related protein,VGR-1,Bmp-6,Vgr1
Accession No.	Uniprot ID:P20722
Uniprot	P20722
GeneID	12161;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Bmp6

Note: This product is for in vitro research use only