

## Bone morphogenetic protein 6

Catalog No: #AP81154



Package Size: #AP81154-1 50ug #AP81154-2 100ug #AP81154-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bone morphogenetic protein 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	BMP-6, VG-1-related protein, VGR-1, Bmp-6, Vgr
Accession No.	Uniprot ID:Q04906
Uniprot	Q04906
GeneID	25644;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Bmp6

Note: This product is for in vitro research use only