

D-beta-hydroxybutyrate dehydrogenase, mitochondrial

Catalog No: #AP81160



Package Size: #AP81160-1 50ug #AP81160-2 100ug #AP81160-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	D-beta-hydroxybutyrate dehydrogenase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	3-hydroxybutyrate dehydrogenase,BDH,Bdh
Accession No.	Uniprot ID:Q80XN0
Uniprot	Q80XN0
GeneID	71911;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Bdh1

Note: This product is for in vitro research use only