

## Transforming growth factor-beta-induced protein ig-h3

Catalog No: #AP81161



Package Size: #AP81161-1 50ug #AP81161-2 100ug #AP81161-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Transforming growth factor-beta-induced protein ig-h3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Beta ig-h3,Kerato-epithelin,RGD-containing collagen-associated protein,RGD-CAP,BIGH3
Accession No.	Uniprot ID:Q15582
Uniprot	Q15582
GeneID	7045;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:TGFBI

Note: This product is for in vitro research use only