

Scavenger receptor cysteine-rich type 1 protein M130

Catalog No: #AP81166



Package Size: #AP81166-1 50ug #AP81166-2 100ug #AP81166-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Scavenger receptor cysteine-rich type 1 protein M130
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	M130,CD163
Accession No.	Uniprot ID:Q2VLG6
Uniprot	Q2VLG6
GeneID	477704;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CD163

Note: This product is for in vitro research use only