

# Beta-1,4-galactosyltransferase 1

Catalog No: #AP81168

Package Size: #AP81168-1 50ug #AP81168-2 100ug #AP81168-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                                                                                                                               |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name       | Beta-1,4-galactosyltransferase 1                                                                                                                              |
| Brief Description  | Recombinant Protein                                                                                                                                           |
| Host Species       | E.coli                                                                                                                                                        |
| Purification       | Greater than 90% by SDS-PAGE                                                                                                                                  |
| Species Reactivity | Ms                                                                                                                                                            |
| Other Names        | Beta-1,4-GalTase 1,UDP-Gal:beta-GlcNAc beta-1,4-galactosyltransferase<br>1,UDP-galactose:beta-N-acetylglucosamine beta-1,4-galactosyltransferase 1,Ggtb,Ggtb2 |
| Accession No.      | Uniprot ID:P15535                                                                                                                                             |
| Uniprot            | P15535                                                                                                                                                        |
| GeneID             | 14595;                                                                                                                                                        |
| Tag Info           | His                                                                                                                                                           |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)                                                                                                        |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:B4galt1

Note: This product is for in vitro research use only