

Carboxypeptidase B2

Catalog No: #AP81171



Package Size: #AP81171-1 50ug #AP81171-2 100ug #AP81171-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Carboxypeptidase B2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Carboxypeptidase U,CPU,Plasma carboxypeptidase B,pCPB,Thrombin-activable fibrinolysis inhibitor,TAFI
Accession No.	Uniprot ID:Q96IY4
Uniprot	Q96IY4
GeneID	1361;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CPB2

Note: This product is for in vitro research use only