

C4b-binding protein alpha chain

Catalog No: #AP81172



Package Size: #AP81172-1 50ug #AP81172-2 100ug #AP81172-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	C4b-binding protein alpha chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	C4bp
Accession No.	Uniprot ID:Q63514
Uniprot	Q63514
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:C4bpa

Note: This product is for in vitro research use only