

## Calpain-1 catalytic subunit

Catalog No: #AP81175



Package Size: #AP81175-1 50ug #AP81175-2 100ug #AP81175-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Calpain-1 catalytic subunit   |
| Brief Description  | Recombinant Protein   |
| Host Species       | E.coli  |
| Purification       | Greater than 90% by SDS-PAGE  |
| Species Reactivity | Bovine  |
| Other Names        | Calcium-activated neutral proteinase 1,CANP 1,Calpain mu-type,Calpain-1 large subunit,Micromolar-calpain,muCANP |
| Accession No.      | Uniprot ID:Q27970   |
| Uniprot            | Q27970  |
| GeneID             | 281661;   |
| Tag Info           | His   |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)  |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CAPN1

Note: This product is for in vitro research use only