

Calpain-5

Catalog No: #AP81177



Package Size: #AP81177-1 50ug #AP81177-2 100ug #AP81177-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Calpain-5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Calpain htra-3,New calpain 3,nCL-3,NCL3
Accession No.	Uniprot ID:O15484
Uniprot	O15484
GeneID	726;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CAPN5

Note: This product is for in vitro research use only