

# B-cell lymphoma 6 protein homolog

Catalog No: #AP81178



Package Size: #AP81178-1 50ug #AP81178-2 100ug #AP81178-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	B-cell lymphoma 6 protein homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	RCJMB04_3d20
Accession No.	Uniprot ID:Q5ZM39
Uniprot	Q5ZM39
GeneID	424912;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: BCL6

Note: This product is for in vitro research use only