

Cadherin-7

Catalog No: #AP81180



Package Size: #AP81180-1 50ug #AP81180-2 100ug #AP81180-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cadherin-7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cad7
Accession No.	Uniprot ID:Q5DWV2
Uniprot	Q5DWV2
GeneID	29162;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cdh7

Note: This product is for in vitro research use only