

## Calpain-1 catalytic subunit

Catalog No: #AP81181



Package Size: #AP81181-1 50ug #AP81181-2 100ug #AP81181-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Calpain-1 catalytic subunit
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Calcium-activated neutral proteinase 1,CANP 1,Calpain mu-type,Calpain-1 large subunit,Micromolar-calpain,muCANP,Canp1,Capa1
Accession No.	Uniprot ID:O35350
Uniprot	O35350
GeneID	12333;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Capn1

Note: This product is for in vitro research use only