

Calpain-2 catalytic subunit

Catalog No: #AP81182



Package Size: #AP81182-1 50ug #AP81182-2 100ug #AP81182-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Calpain-2 catalytic subunit
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Calcium-activated neutral proteinase 2,CANP 2,Calpain M-type,Calpain-2 large subunit,Millimolar-calpain,M-calpain
Accession No.	Uniprot ID:Q92178
Uniprot	Q92178
GeneID	395963;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CAPN2

Note: This product is for in vitro research use only