

Transcription factor BTF3

Catalog No: #AP81186



Package Size: #AP81186-1 50ug #AP81186-2 100ug #AP81186-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transcription factor BTF3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Nascent polypeptide-associated complex subunit beta,NAC-beta,RNA polymerase B transcription factor 3
Accession No.	Uniprot ID:Q64152
Uniprot	Q64152
GeneID	218490;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Btf3

Note: This product is for in vitro research use only