

Cadherin-11

Catalog No: #AP81187



Package Size: #AP81187-1 50ug #AP81187-2 100ug #AP81187-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cadherin-11
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	OSF-4,Osteoblast cadherin,OB-cadherin,Cad-11
Accession No.	Uniprot ID:P55288
Uniprot	P55288
GeneID	12552;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cdh11

Note: This product is for in vitro research use only