

## Carbonic anhydrase 1

Catalog No: #AP81188



Package Size: #AP81188-1 50ug #AP81188-2 100ug #AP81188-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Carbonic anhydrase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Carbonate dehydratase I,Carbonic anhydrase I,CA-I
Accession No.	Uniprot ID:Q1LZA1
Uniprot	Q1LZA1
GeneID	510801;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CA1

Note: This product is for in vitro research use only