

Calpain-3

Catalog No: #AP81195

Package Size: #AP81195-1 50ug #AP81195-2 100ug #AP81195-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Calpain-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Calcium-activated neutral proteinase 3,CANP 3,Calpain L3,Calpain p94,Muscle-specific calcium-activated neutral protease 3
Accession No.	Uniprot ID:Q64691
Uniprot	Q64691
GeneID	12335;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Capn3

Note: This product is for in vitro research use only