

Macrophage-capping protein

Catalog No: #AP81196



Package Size: #AP81196-1 50ug #AP81196-2 100ug #AP81196-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Macrophage-capping protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Actin regulatory protein CAP-G
Accession No.	Uniprot ID:Q6AYC4
Uniprot	Q6AYC4
GeneID	297339;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Capg

Note: This product is for in vitro research use only