

Caspase-14

Catalog No: #AP81198

Package Size: #AP81198-1 50ug #AP81198-2 100ug #AP81198-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Caspase-14
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CASP-14,Mini-ICE,MICE
Accession No.	Uniprot ID:O89094
Uniprot	O89094
GeneID	12365;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Casp14

Note: This product is for in vitro research use only