

C5a anaphylatoxin chemotactic receptor C5L2

Catalog No: #AP81203



Package Size: #AP81203-1 50ug #AP81203-2 100ug #AP81203-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C5a anaphylatoxin chemotactic receptor C5L2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	G-protein coupled receptor 77,C5I2
Accession No.	Uniprot ID:Q695P6
Uniprot	Q695P6
GeneID	445269;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gpr77

Note: This product is for in vitro research use only