

Bone morphogenetic protein 1

Catalog No: #AP81207



Package Size: #AP81207-1 50ug #AP81207-2 100ug #AP81207-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Bone morphogenetic protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	BMP-1,Mammalian tolloid protein,mTld,Procollagen C-proteinase,PCP,PCOLC
Accession No.	Uniprot ID:P13497
Uniprot	P13497
GeneID	649;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BMP1

Note: This product is for in vitro research use only