

Calpain-6

Catalog No: #AP81210



Package Size: #AP81210-1 50ug #AP81210-2 100ug #AP81210-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Calpain-6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Capa6
Accession No.	Uniprot ID:O35646
Uniprot	O35646
GeneID	12338;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Capn6

Note: This product is for in vitro research use only