

Baculoviral IAP repeat-containing protein 2

Catalog No: #AP81212

Package Size: #AP81212-1 50ug #AP81212-2 100ug #AP81212-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Baculoviral IAP repeat-containing protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cellular inhibitor of apoptosis 1,C-IAP1,Inhibitor of apoptosis protein 2,mlAP2,RING-type E3 ubiquitin transferase BIRC2
Accession No.	Uniprot ID:Q62210
Uniprot	Q62210
GeneID	11797;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Birc2

Note: This product is for in vitro research use only