

Caspase-1

Catalog No: #AP81213

Package Size: #AP81213-1 50ug #AP81213-2 100ug #AP81213-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Caspase-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CASP-1,Interleukin-1 beta convertase,IL-1BC,Interleukin-1 beta-converting enzyme,ICE,p45,Il1bc
Accession No.	Uniprot ID:P29452
Uniprot	P29452
GeneID	12362;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Casp1

Note: This product is for in vitro research use only