

Coiled-coil domain-containing protein 104

Catalog No: #AP81214



Package Size: #AP81214-1 50ug #AP81214-2 100ug #AP81214-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Coiled-coil domain-containing protein 104
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Accession No.	Uniprot ID:Q4V8E4
Uniprot	Q4V8E4
GeneID	289859;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ccdc104

Note: This product is for in vitro research use only