

## Carbonic anhydrase 4

Catalog No: #AP81217



Package Size: #AP81217-1 50ug #AP81217-2 100ug #AP81217-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Carbonic anhydrase 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Carbonate dehydratase IV,Carbonic anhydrase IV,CA-IV
Accession No.	Uniprot ID:Q95323
Uniprot	Q95323
GeneID	280741;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CA4

Note: This product is for in vitro research use only