

Calcitonin gene-related peptide 2

Catalog No: #AP81218



Package Size: #AP81218-1 50ug #AP81218-2 100ug #AP81218-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Calcitonin gene-related peptide 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Beta-type CGRP, Calcitonin gene-related peptide II, CGRP-II
Accession No.	Uniprot ID:P10093
Uniprot	P10093
GeneID	171519;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Calc

Note: This product is for in vitro research use only